



CERTIFICATION

AOAC Research Institute *Performance Tested Methods*SM

Certificate No.
100903

The AOAC Research Institute hereby certifies the method known as:

VitaFast® Folic Acid Microbiological Microtiter Plate Test for the Determination of Folic Acid

manufactured by
ifp Institut für Produktqualität GmbH
Wagner-Régeny-Str. 8
12489 Berlin
Germany

distributed by
R-Biopharm AG
An der neuen Bergstraße 17
64297 Darmstadt
Germany

This method has been evaluated in the AOAC Research Institute *Performance Tested Methods*SM Program and found to perform as stated in the applicability of the method. This certificate indicates an AOAC Research Institute Certification Mark License Agreement has been executed which authorizes the manufacturer to display the AOAC Research Institute *Performance Tested Methods*SM certification mark on the above-mentioned method for the period below. Renewal may be granted by the Expiration Date under the rules stated in the licensing agreement.

A handwritten signature in black ink that reads "Scott Coates".

Scott Coates, Senior Director
Signature for AOAC Research Institute

Issue Date	December 2, 2022
Expiration Date	December 31, 2023

AUTHORS ORIGINAL VALIDATION: Jessica Kerr and Kurt Johnson MODIFICATION FEBRUARY 2017: Wolfgang Weber	SUBMITTING COMPANY R-Biopharm Inc. 7950 Old US 27 South Marshall, MI 49068	Current Sponsor R-Biopharm AG An der neuen Bergstraße 17 64297 Darmstadt Germany
--	--	---

METHOD NAME VitaFast® Vitamin Folic Acid Microbiological Microtiter Plate Test for the Determination of Folic Acid	CATALOG NUMBER P1001
--	--------------------------------

INDEPENDENT LABORATORY University of Guelph Laboratory Services Division 95 Stone Road West Guelph, ON N1G 2Z4 Canada	AOAC EXPERTS AND PEER REVIEWERS Sneh Bhandari ^{1,3} , Michael Rychlik ² ¹ Silliker Laboratories, Illinois, USA ² Technische Universität München, GERMANY ³ Modifications: February 2017
--	---

APPLICABILITY OF METHOD Target analyte – B vitamin folic acid Matrixes – (1 g) - Cereals, infant formula, processed meats, multivitamin pills, powders, beverages like fruit juice & milk.	Performance claims – The performance characteristics of VitaFast® Folic acid meet the following specifications: 1) Time required for completion of the sample extraction was 2 hours and less than 48 hours for the test implementation. 2) The test kit components are stable as indicated on the test kit labels (shelf life is 12 months). 3) Analytical Sensitivity was found at LOD 0.018 µg / 100 g (ml) as measured by 190 blank samples from 10 different lots. LOQ was set at 0.16 µg Folic acid / 100 g sample, which corresponds to standard 1 of the curve. 4) Accuracy was investigated by analysis of reference materials from five proficiency programs, and by commercial product analysis and spike recovery studies. In general recovery was within acceptable limits. 5) The VitaFast test kit was shown to have a high degree of precision, with inter-assay variances below 10 % for all matrixes. 6) The VitaFast plate test is not sensitive to temperature changes between 36 and 38 °C, incubation time between 44 and 52 hours, or assay medium volumes between 145 and 155 µl.
---	--

ORIGINAL CERTIFICATION DATE October 28, 2009	CERTIFICATION RENEWAL RECORD Renewed annually through December 2023.
--	--

METHOD MODIFICATION RECORD 1. February 2017 Level 2	SUMMARY OF MODIFICATION 1. Manufacturing location change.
---	---

Under this AOAC <i>Performance Tested Methods</i> SM License Number, 100903 this method is distributed by: R-Biopharm AG	Under this AOAC <i>Performance Tested Methods</i> SM License Number, 100903 this method is distributed as: VitaFast® Vitamin B2 (Riboflavin) Microbiological Microtiter Plate Test for the Determination of Vitamin B2 Riboflavin
--	---

PRINCIPLE OF THE METHOD (1)

Folic acid is extracted from the sample and the extract is diluted. The diluted extract and the folic acid assay - medium are pipetted into the wells of a microtiter plate which are coated with Lactobacillus rhamnosus. The growth of L. rhamnosus is dependent on the supply of folic acid. Following the addition of folic acid as a standard or as a compound of the sample, the bacteria grow until the vitamin is consumed. The incubation is done in the dark at 37 °C (98.6 °F) for 44 - 48 h. The intensity of metabolism or growth in relation to the extracted folic acid is measured as turbidity and compared to a standard curve. The measurement is done using a microtiter plate reader at 610 - 630 nm (alternatively at 540 - 550 nm).

DISCUSSION OF THE VALIDATION STUDY (1)

The VitaFast Folic acid test kit is calibrated according to a standard curve of five standard concentrations, using 4-parameter fitting software. The curve shown in figure 1 is typical. Variation within the curve is consistently minor, at a level of variance below 10 %. Stability is also demonstrated over the entire shelf life of the product, and regular quality tests ensure this is true for all lots produced.

In the independent laboratory study, VitaFast produced a mean result for the NIST reference material approximately 7% higher than the upper limit of uncertainty, whereas the AOAC reference method 944.12 produced a result of almost 100% the declared value. Although the VitaFast method did produce a result slightly above the specification of the sample, still demonstrates reasonable accuracy and excellent precision, as shown by the RSD.

Furthermore, VitaFast produced a result for the cereal sample very close to the label claim (104% recovery). AOAC method 944.12 also agreed closely with the label claim for the cereal product, yielding a result of 93% recovery.

Table 7 Intra-assay variance of food samples (triplicate analysis per sample dilution) (1)				
Sample description conc. indicated on label in µg / 100 g (ml)	Dilution factor	Mean result in µg / 100 g (ml)	Mean result of dilutions in µg / 100 g (ml)	Coefficient of variation in %
Cereal flakes with fruits (261)	2000	254.6	275.7	7.6
	1500	296.6		
	1000	276.0		
Diet potato soup (100)	600	117.4	123.2	6.5
	400	132.4		
	200	119.9		
Tomato soup (200)	1200	108.0	111.8	9.6
	800	103.4		
	400	123.9		
Infant formula (15)	100	29.1	28.2	4.1
	60	28.6		
	40	26.9		
Semolina with milk (62)	100	101	104	6.9
	150	96		
	200	104		
	300	113		
Hospital milk drink (60)	100	81	80	1.8
	200	79		
Baby formula (140)	200	91	90	1.6
	400	89		

Table 8. Intra-assay variance of food samples (1)			
Sample description	Conc. indicated on label in µg / 100 g	Mean result of dilutions in µg / 100 g	Coefficient of variation in %
Dextrose powder RM – Vit001 internal Reference material	200	204 (n=4)	6.8
		211 (n=3)	8.4
		205 (n=4)	6.2
		204 (n=4)	6.6
Ham sausage	200	415 (n=3)	2.1
Cake	150	119 (n=4)	6.0
Cereal (tested by four different groups – internal training)	500	617 (n=6)	7.5
		557 (n=6)	
		665 (n=6)	
		642 (n=6)	
Pineapple tangerine drink	30	28 (n=2)	1.7
		29 (n=6)	
Multivitamin juice	30	40 (n=2)	1.2
		41 (n=6)	
Multivitamin nectar	30	33 (n=2)	2.8
		34 (n=2)	
		35 (n=4)	

REFERENCES CITED

1. Kerr, Jessica, and Johnson, Kurt., Evaluation of the VitaFast® Folic Acid: Microbiological Microtiter Plate Test for the Determination of Folic Acid, AOAC *Performance Tested Methods*SM certification number 100903.
2. AOAC Research Institute Validation Outline for VitaFast® Folic Acid: Microbiological Microtiter Plate Test for the Determination of Folic Acid, Approved – October 2009.
3. Higdon, J. Linus Pauling Institute, Oregon State University. 2008. Micronutrient Information Center: Vitamin B12. [Accessed online January 25, 2008: <http://lpi.oregonstate.edu/infocenter/vitamins/pa/>]
4. Food and Nutrition Board, Institute of Medicine. 1998. Dietary Reference Intakes for Thiamin, Riboflavin, Niacin, Vitamin B6, Folate, Vitamin B12, Pantothenic Acid, Biotin and Choline. [Accessed online April 20, 2009: http://www.nap.edu/catalog.php?record_id=6015#toc]
5. European Food Information Council. 2006. Vitamins: What they do and where to Find them. [Accessed online april 15, 2009: <http://www.eufic.org/article/en/page/MARCHIVE/expid/miniguide-vitamins/#9>]
6. Health Canada. 2005. Addition of Vitamins and Minerals to Foods. [Accessed online January 25, 2008: http://www.hc-sc.gc.ca/fn-an/nutrition/vitamin/fortification_final_doc_1-eng.php]
7. U.S. Food and Drug Administration. 2009. Fortify Your Knowledge about Vitamins. [Accessed online April 15, 2009: <http://www.fda.gov/consumer/updates/vitamins111907.html>]
8. Ball, G.F.M. 1998. Bioavailability and Analysis of Vitamins in Foods. Chapman & Hall.
9. Allen, L.H. 2003. B Vitamins: Proposed Fortification Levels for Complementary Foods for Young Children. J. Nutr. 133: 3000S-3007S
10. Weber, W., Modification Report for Location Change, AOAC *Performance Tested Methods*SM certification number 100903. Approved February 2017.